



Biomedical implications of H₂ enriched drinking water

Unisense hydrogen sensors used for detection of H₂ in biomedical research and in production control of H₂ enriched water

The application note is based on the research and articles by:
Kamimura et al. (2011)

The application note is written by:
Daugaard et al., Unisense

H₂ sensor characteristics:

Sensor for dissolved and gaseous hydrogen

Modifications and customizations available

Detection range <50 nM to saturated H₂ water

Response time in seconds

Real-time data

Easy calibration and linear response

H₂ enriched drinking water

Drinking water enriched with dissolved hydrogen has been reported to have a positive effect on many health and well being conditions. This has made a market, predominately in Japan, to commercialize H₂-enriched water for health and body recovery processes.

Unisense hydrogen microsensors have been used to document the possible positive effect of H₂-enriched drinking water as a biomedical treatment. In the study by Kamimura et al. (2011) they found that obese rats that received H₂ water were less affected by oxidative stress in connection with obesity compared to rats that received normal water.

Using a needle type H₂ microsensor, the H₂ concentration was monitored in rat liver. The measurements showed that H₂ was accumulated in liver of fed rats but not in fasted rats (Kamimura personal communication).




Figure 1: Hydrogen concentration measured in rat liver of a fed or overnight-fasted rat liver. Arrow indicates when the rats received H₂ enriched drinking water.



The study also showed that a fatty liver from the rats drinking H₂ rich water contained less oxidative stress marker proteins and looked healthier compared to the fat rats, that did not receive H₂ enriched water. Hydrogen's anti-oxidative effect on the liver is speculated to explain the difference.

Unisense H₂ microsensors are used by manufacturers of commercial hydrogen enriched water to control and monitor the dissolved H₂ concentration.



Examples of commercially available hydrogen enriched drinking water




Figure 1: Hydrogen concentration in open tubes with water and glycogen respectively after bubbling with H₂. Data are based on three replicates, * P<0.05, *** P<0.001.

Suggested products

- Various customizations of the H₂ sensor



Version: June 2024

Related publications

Kamimura et al. (2011), Molecular Hydrogen Improves Obesity and Diabetes by Inducing Hepatic FG21 and Stimulating Energy Metabolism in db/db Mice, Obesity, Vol. 19: 1396-1403.

Terasaki et al. (2012) Hydrogen therapy attenuates irradiation-induced lung damage by reducing oxidative stress. Am J Physiol Lung Cell Mol Physiol, Vol. 301: L415-L426.

Sun et al. (2011) The protective role of hydrogen-rich saline in experimental liver injury in mice. Journal of Hepatology, Vol. 54: 471-480.

Yamamoto et al. (2019) Hydrogen gas distribution in organs after inhalation: Real-time monitoring of tissue hydrogen concentration in rat. Scientific Reports, Vol 9, Issue 1. 1255.

Terasaki et al. (2019) Molecular hydrogen attenuates gefitinib-induced exacerbation of naphthalene-evoked acute lung injury through a reduction in oxidative stress and inflammation, Laboratory Investigation.

Iketani et al. (2018) Administration of hydrogen-rich water prevents vascular aging of the aorta in LDL receptor-deficient mice, Scientific Reports, Vol 8, Issue 1, p. 1-11.

He et al. (2017) Image-Guided Hydrogen Gas Delivery for Protection from Myocardial Ischemia-Reperfusion Injury via Microbubbles, ACS Applied Materials and Interfaces, Vol 9, Issue 25, p. 21190-21199.