



Unisense microsensors in water splitting studies

Hydrogen and oxygen microsensors for
direct and real-time detection of dissolved
gas in photocatalytic/electrochemical studies

The application note is based on
the research and articles by:

Sick et al. (2018)

Wang et al. (2016)

The application note is written by:

Daugaard et al., Unisense

H₂ and O₂ sensor characteristics:

Direct gas and liquid phase
measurements

Real-time data
- no need for GC or MS

Detection range
<50 nM to saturation

Response time in seconds

Customizations available

Simple two-point calibration

Selected microsensor applications:

Hydrogen evolved from light-driven water splitting at the surface of a semiconductor material

Jonathan Kampmann and coworkers from the University of Munich have used a Unisense H₂-NP low range microsensor to detect the nanomolar concentrations of hydrogen evolved from light-driven water splitting at the surface of a semiconductor material. They used a covalent organic framework (COF) as photoelectrode and the material was held at a potential of 0.4 V versus a reversible hydrogen electrode. During light illumination, the generated charge reduced water in the cuvette and hydrogen was produced and detected continuously with the H₂ microsensor (Figure 1).

You can read more in the article by Sick et al. (2018) Oriented Films of Conjugated 2D Covalent Organic Frameworks as Photocathodes for Water Splitting. *J. Am. Chem. Soc.* 2018, 140, 2085–2092.




Figure 1: Microsensor measurement of the H₂ produced when the sample was illuminated with sunlight (left). H₂-NPLR microsensor in the cuvette containing the active semiconducting material (orange substrate) and water (right).

Data and photo kindly provided by Jonathan Kampmann.



O₂

N₂O

H₂S

NO

H₂

pH

Redox

Temp

EP

Hydrogen production upon illumination of a polymer dot catalyst

Dr. Haining Tian and his research group from Uppsala University have been investigating the potential of organic semiconducting polymer dots (Pdots) to serve as a photocatalyst in light-driven hydrogen production from water splitting. In the experiment the researchers used a H₂-NP low range sensor with a detection limit of 50 nM to pierce through the lid of a cuvette and then measure the hydrogen evolution in an aqueous solution containing Pdots (Figure 2). The researchers showed that the Pdots worked very well as photocatalysts for light-driven water splitting.






Figure 1: H₂ generation over time measured with the Unisense H₂-NPLR microsensor (left) and experimental setup showing the hydrogen microsensor pierced through the lid into the sample measuring the H₂ generation (above).

Data and photo kindly provided by Dr. Haining Tian.

For further reading please see the article: Wang et al. (2016) Organic Polymer Dots as Photocatalysts for Visible Light-Driven Hydrogen Generation. *Angew. Chem. Int. Ed.* 2016.

Suggested products

- Sensors
- Single or MultiChannel UniAmp
- SensorTrace Rate



Related publications

Salomao et al. (2019) Photoelectrochemical hydrogen production from water splitting using heterostructured nanowire arrays of Bi₂O₃/BiAl oxides as photocathode. *Solar Energy Materials and Solar Cells* 2019.

Alqathani et al. (2019) Heteroepitaxy of GaP on silicon for efficient and cost-effective photoelectrochemical water splitting. *J. Mater. Chem. A*, 7, 8550-8558.

Tian et al (2018) Hydrogen evolution by a photoelectrochemical cell based on a Cu₂O-ZnO-[FeFe] hydrogenase electrode. *J. Photochem. Photobiol. A Chem.*, Vol 366, 27-33.

Sherman et al. (2016) A Dye-Sensitized Photoelectrochemical Tandem Cell for Light Driven Hydrogen Production from Water, *J. Am. Chem. Soc.* 138, 16745-16753.