

## OXY-METER AND pH/mV-METER USER MANUAL



# **OXY-METER AND pH/mV-METER USER MANUAL**

Copyright © 2018 · Unisense A/S

Version July 2018

# **OXY-METER AND pH/mV-METER USER MANUAL**

UNISENSE A/S



# TABLE OF CONTENTS

CONGRATULATIONS WITH YOUR NEW PRODUCT!.....	6
SUPPORT, ORDERING, AND CONTACT INFORMATION	6
OVERVIEW .....	8
SPECIFICATIONS .....	10
GETTING STARTED .....	11
ELECTRODES	11
OPERATING ENVIRONMENT	11
WARRANTY AND LIABILITY.....	12
REPAIR OR ADJUSTMENT	12

# CONGRATULATIONS WITH YOUR NEW PRODUCT!

## **SUPPORT, ORDERING, AND CONTACT INFORMATION**

If you wish to order additional products or if you encounter any problems and need scientific/technical assistance, please do not hesitate to contact our sales and support team. We will respond to your inquiry within one working day.

E-mail: [sales@unisense.com](mailto:sales@unisense.com)

Unisense A/S  
Tueager 1  
DK-8200 Aarhus N, Denmark  
Tel: +45 8944 9500  
Fax: +45 8944 9549

Further documentation and support is available at our website  
[www.unisense.com](http://www.unisense.com).



# OVERVIEW

The Unisense one-channel amplifiers are dedicated picoammeters and high-impedance millivoltmeters with a built-in A/D-converter, designed for use with Unisense Clark-type micro and mini sized sensors and electrodes.

Unisense make dedicated standard one-channel amplifiers for all our sensors and electrodes. We also make custom-built amplifiers where the customer selects parameters such as channel composition and focus sensitivity



*Front and back views of OXY-Meter (top) and pH/mV-Meter (bottom)*



# SPECIFICATIONS

<b>INPUT CONNECTORS</b>	
Sensor	Lemo
Reference (mV channels)	Banana plug
<b>AMPLIFIER</b>	
Types	O <sub>2</sub> , pH/redox
Number	1
Polarization	Preset
Input range	± 4500 pA, ± 4500 mV, depending on type, or user-defined
Input impedance (mV )	> 10 <sup>13</sup> ohm
<b>RESPONSE AND CONTROL</b>	
Response time (90%)	200 ms
<b>POWER</b>	
Power supply	USB
<b>A/D-CONVERTER</b>	
Resolution	12 bit
Sampling frequency	Up to 1.000 Hz
<b>INTERFACE</b>	
Digital	USB
<b>DATA ACQUISITION SOFTWARE</b>	
SensorTrace Logger	Included
Other Unisense software	Optional: SensorTrace Profiling
Data rate (max)	10 samples/sec
Software data output	CSV formatted file (MS Excel compatible)
<b>PHYSICAL</b>	
Dimensions ((W x D x H)	113 x 125 x 50 mm
Weight	720 g
<b>ENVIRONMENTAL</b>	
Temperature range	-10 to +70°C

# GETTING STARTED

1. Install the software
2. Connect the amplifier to the USB port (USB 2.0) of the computer
3. Connect the sensor/electrodes. If relevant, ground the set-up by connecting the provided blue grounding cable
4. Start the software (consult the software manual for directions)
5. The sensor signal(s) should now be visible on the computer screen and calibrations and measurements can be performed.

## ELECTRODES

The Unisense dedicated amplifiers are designed for use with Unisense microsensors and electrodes with LEMO connectors. If you have a microsensor or electrode with BNC connector, Unisense can supply a BNC to LEMO adaptor cable.

### **IMPORTANT:**

*Please consult the separate sensor manual for instructions about the oxygen sensor, prepolarization, and calibration.*

## OPERATING ENVIRONMENT

The one-channel amplifiers are designed for laboratory use. Care must be taken to protect the amplifiers in humid, dusty, and corrosive environments (e.g. field conditions) since it may affect the performance, cause damage, or shorten lifetime. If the amplifier is used outdoors, it should be protected in a plastic bag or similar.

# WARRANTY AND LIABILITY

Microsensors are consumable items. The sensors are tested thoroughly before packaging and shipment.

Warranty for sensors:

N<sub>2</sub>O sensors: 60 days from shipment.

Standard Oxygen sensors: 180 days from shipment.

All other sensors excluding special sensors: 90 days from shipment.

If, within the above specified period, the sensor(s) fail to perform according to the specifications, Unisense will replace the sensor(s) free of charge.

Unisense will only replace dysfunctional sensors if they have been tested according with the instructions in the manual upon receipt of the sensor(s). The warranty does not include repair or replacement necessitated by accident, neglect, misuse, unauthorized repair, or modification of the product. Physical damage to the tip of the sensor is not covered by the warranty.

## Liability

In no event will Unisense be liable for any direct, indirect, consequential or incidental damages, including lost profits, or for any claim by any third party, arising out of the use, the results of use, or the inability to use this product.

## REPAIR OR ADJUSTMENT

Sensors and electrodes cannot be repaired. Other equipment that is not covered by the warranty will, if possible, be repaired by Unisense with appropriate charges paid by the customer. In case of return of equipment please contact us for a return authorization.

For further information please see the documents "Conditions for Sale and Delivery for Unisense" and "Warranty and Shipping Information".



**UNISENSE**

UNISENSE, DENMARK

[www.unisense.com](http://www.unisense.com) · [info@unisense.com](mailto:info@unisense.com)