

## Temperature response of Unisense oxygen microsensors

The signal of an oxygen sensor is temperature dependent, as factors influencing the sensor signal varies with temperature. Within the temperature range from 5-30 C the signal  $S$  of a Unisense oxygen sensor, at the temperature  $T$  and partial pressure  $P$  can be described as (1):

$$S(P, T) = Z + a(T)P \quad (1)$$

The zero current  $Z$  for a Unisense sensor is approximately temperature independent. The coefficient  $a$  is temperature dependent as described by equation 2:

$$a(T) = Ae^{kT} \quad (2)$$

By substituting 2 into 1 we get the following expression

$$S(P, T) = Z + Ae^{kT} P \quad (3)$$

The sensor signal at a known temperature  $T_0$  and partial pressure  $P_0$  is

$$S(P_0, T_0) = Z + Ae^{kT_0} P_0 \quad (4)$$

Which can be rearranged to

$$A = \frac{S(P_0, T_0) - Z}{e^{kT_0} P_0} \quad (5)$$

Substituting 5 into 3 and 4, we can rewrite these to

$$S(P, T) = Z + (S(P_0, T_0) - Z) \frac{P}{P_0} e^{k(T-T_0)} \quad (6)$$

or

$$\ln(S(P, T) - Z) = \ln(S(P_0, T_0) - Z) + \ln \frac{P}{P_0} + k(T - T_0) \quad (7)$$

By plotting the signal  $S(P, T)$  as a function of temperature  $T$  a fixed partial pressure  $P$  in equation 7,  $k$  can be determined as the slope to a linear fit of the plot.

$$P = \frac{S(P, T) - Z}{S(P_0, T_0) - Z} P_0 e^{k(T_0 - T)} \quad (8)$$

Once  $k$  is determined it is from (8) possible to determine the partial pressure  $P$  from sensor signals  $S$  and temperature measurements  $T$ , providing the zero current  $Z$ , and a standard signal at a known temperature and partial pressure  $S(P_0, T_0)$ , is known.

The temperature coefficient  $k$  can vary from sensor to sensor due to differences in internal dimensions and membrane thickness, although it for most sensors is between  $0.015 \text{ T}^{-1}$  and  $0.025 \text{ T}^{-1}$ . The temperature coefficient has to be determined individually for each sensor if an exact correction is needed.