




UNISENSE



The Microsensor Company



Demonstration of Activity Calculation - Software




February 2022
Tage Dalsgaard


A few rules before we get started



1. Please turn off your microphone



2. Questions: During lecture please use chat.
After lecture you can unmute and ask.



Very application-specific questions may be better answered in a private session afterwards.

You will get access to all the presentations as PDF's + recordings shortly after the workshop.

Online Biogeochemistry Workshop



Wednesday 23 February

14:30-15:30 CET - Introduction to Microsensors

15:45-16:45 CET - Lab-based Studies

Thursday 24 February

14:30-15:30 CET - Field Studies

15:45-16:30 CET - Demonstration of Field Microprofiling System

16:45-17:30 CET - Demonstration of Activity Calculation - Software

Profile analysis

Calculation of production and consumption rates from microprofiles

Based on the method from:

Berg, P., N. Risgaard-Petersen, and S. Rysgaard. 1998. Interpretation of measured concentration profiles in sediment pore water. *Limnol. Oceanogr.* 43: 1500–1510.

Limnol. Oceanogr., 43(7), 1998, 1500–1510
© 1998, by the American Society of Limnology and Oceanography, Inc.

Interpretation of measured concentration profiles in sediment pore water

*Peter Berg*¹

Department of Environmental Sciences, University of Virginia, Charlottesville, Virginia 22903

Nils Risgaard-Petersen


Institute of Biological Science, Department of Microbial Ecology, University of Aarhus, Ny Munkegade, Bldg. 540, DK-8000 Aarhus C, Denmark; and National Environmental Research Institute, Department of Lake and Estuarine Ecology, Vejlsovej 25, DK-8600 Silkeborg, Denmark

Søren Rysgaard

National Environmental Research Institute, Department of Lake and Estuarine Ecology, Vejlsovej 25, DK-8600 Silkeborg, Denmark

Abstract

A robust numerical procedure for biogeochemical interpretation and analysis of measured concentration profiles of solutes in sediment pore water has been developed. Assuming that the concentration-depth profile represents a




UNISENSE


Demo profiles


Demonstration profiles installed with *SensorTrace Suite*

- ProfilingDemo.upro – Read only
 - \Documents\Unisense Data\Demo Experiments
- Profiling application note
 - Unisense.com



The screenshot shows a web browser window displaying the Unisense website. The URL is unisense.com/applications/sediment-profiling-analysis/. The navigation menu includes 'Products', 'Knowledge' (circled in red), 'Support', 'About us', and 'Order now'. Below the navigation, there is a breadcrumb trail: 'Home > Application Notes > Sediment Profiling Analysis'. The main content area features a large image of a sediment profiler in a shallow, brownish water body. Overlaid on the image is the text: 'CALCULATIONS OF OXYGEN CONSUMPTION RATE USING SENSORTRACE PROFILING' and 'Sediment Profiling Analysis'. Below the image, there is a table listing application notes.

Name	Filetype
 Profiling Application Note (PDF)	PDF, 2 MB



The cover page features the Unisense logo at the top right. Below it is a large image of a sediment profiler in a shallow, brownish water body. In the bottom left corner, there is a small icon of a graph and a gear. The main title is 'Calculations of oxygen consumption rate using SensorTrace Profiling'. Below the title, there are two columns of text: 'Abstract' and 'Material and Method'.

Abstract
This note will explain how to use SensorTrace Profiling to quantify the consumption rate of oxygen as well as the oxygen exchange rate across

Material and Method
Sediment cores were collected by hand at Aggersund, Limfjorden, Denmark in August 2015 at a water depth of about 30 cm and brought to the


Profile analysis

Calculation of production and consumption rates from microprofiles

Select profile, depth range and boundary conditions

Specify diffusion coefficient and porosity

Output with rates and statistic parameters



O₂ production



O₂ consumption

Fitted profile

Microprofiles - Theory

Example:

Biofilm with O_2 production in top layer, O_2 consumption below




Qualitative information from profiles

1. Straight line: No net consumption or production – only diffusional transport.
2. Concave: Net production (e.g. O_2 production by photosynthesis, H_2S production by sulphate reduction)
3. Convex: Net consumption (e.g. Respiration, oxidation of reduced compounds)

Microprofiles - Theory

1-dimensional system – Steady state



Quantitative Information from microprofiles

1. **C(x) - Concentration:** Penetration depth, overlapping zones
2. **dC(x)/dx – Flux:** Into sediment, within sediment


D_w = Diffusion coefficient in water
 $\phi(x)$ = Porosity
 $D_s(x)$ = Diffusivity
 $D_s(x) = D_w \times \phi(x)$ (simplest form)

Diffusion in pure water and sediment

$$J_w = -D_w \frac{dC(x)}{dx}$$


$$J_{sed} = -\phi(x)D_s(x) \frac{dC(x)}{dx}$$

Pure water




100% space for diffusion
Random movement

Porosity




Example:
Porosity = 0.8 =>
80% space for diffusion


Tortuosity



Longer diffusion path
 $D_s(x) = D_w \times \text{Porosity}$




Sediment diffusion coefficient equation



Microprofiles - Theory

1-dimensional system – Steady state




Quantitative Information from microprofiles

- C(x) - Concentration:** Penetration depth, overlapping zones
- dC(x)/dx – Flux:** Into sediment, within sediment
- d²C(x)/dx² – Production and consumption:** Activity distribution within the sediment

- D_w = Diffusion coefficient in water
- $\phi(x)$ = Porosity
- $D_s(x)$ = Diffusivity
- $D_s(x) = D_w \times \phi(x)$ (simplest form)

Microprofiles from different instruments

Lab Microprofiling System




Lab conditions

Field Microprofiling System



Down to 10 m water depth

In situ Microprofiling System





Down to 6000 m water depth



SensorTrace Profiling will record or import profiles from all system for analysis



Find the sediment surface

Stereo microscope while measuring



Find the sediment surface

Locate surface after measuring



Setting depth interval and boundary conditions

Settings and Analysis

Sensor 1 - O₂

Profile 2.1

Analysis 6

Boundary Conditions

Bottom conc + bottom flux

Bottom conc: 0.27

Bottom flux: 0.00

Intervals and zones

Start depth(μm): 0

Max zones: 7

End depth(μm): 1200

Min width (μm): 100

Production


Consumption

Observed

Simulated

DS & Theta

Analyze



DS and Theta

Settings

Zones with different diffusion parameters: 1

Calculate

D0 coefficient free water (10⁻⁵ cm² s⁻¹): O2 Table 2.104

Sediment diffusion coefficient equation: DS = Por*D

Start	End	Porosity	DS
-1500	1300	1	2.104E-05

Confirm Cancel

Enter DS manually

Diffusion coefficient

- Enter the diffusion coefficient for O₂ in water
- Depends on salinity and temperature
- As many depth intervals as you like

DS and Theta ✕

Settings

Zones with different diffusion parameters:


Calculate

D0 coefficient: free water(10⁻⁵ cm² s⁻¹)

Sediment diffusion coefficient equation: DS = Por*D

Start	End	Porosity	DS
-1500	0	1	2.038E-05
0	1200	0.8	1.63E-05

Confirm
Cancel
Enter DS manually



DATA-TABLE 1 by Niels Ramsing & Jens Gundersen

Diffusion coefficient for oxygen at different temperatures and salinities of seawater

Units: 10⁻⁵ cm² s⁻¹

Salinity (‰)	Temperature (°C)																
	0.0	1.0	2.0	3.0	4.0	5.0	6.0	7.0	8.0	9.0	10.0	11.0	12.0	13.0	14.0	15.0	16.0
0.0	1.1041	1.1465	1.1899	1.2344	1.2798	1.3261	1.3734	1.4214	1.4702	1.5198	1.5700	1.6209	1.6723	1.7243	1.7769	1.8300	1.8836
1.0	1.1026	1.1448	1.1881	1.2324	1.2777	1.3239	1.3709	1.4188	1.4675	1.5169	1.5669	1.6176	1.6689	1.7208	1.7732	1.8261	1.8796
2.0	1.1011	1.1432	1.1863	1.2305	1.2756	1.3216	1.3685	1.4162	1.4647	1.5140	1.5639	1.6144	1.6656	1.7173	1.7695	1.8223	1.8756
3.0	1.0996	1.1415	1.1845	1.2285	1.2735	1.3193	1.3661	1.4137	1.4620	1.5111	1.5608	1.6112	1.6622	1.7137	1.7658	1.8185	1.8717
4.0	1.0981	1.1399	1.1827	1.2266	1.2714	1.3171	1.3637	1.4111	1.4593	1.5082	1.5578	1.6080	1.6588	1.7102	1.7622	1.8147	1.8677
5.0	1.0966	1.1382	1.1809	1.2246	1.2693	1.3148	1.3613	1.4085	1.4566	1.5053	1.5547	1.6048	1.6555	1.7068	1.7586	1.8109	1.8638
6.0	1.0950	1.1366	1.1792	1.2227	1.2672	1.3126	1.3589	1.4060	1.4539	1.5024	1.5517	1.6017	1.6522	1.7033	1.7550	1.8072	1.8599
7.0	1.0935	1.1350	1.1774	1.2208	1.2651	1.3104	1.3565	1.4034	1.4512	1.4996	1.5487	1.5985	1.6489	1.6998	1.7514	1.8034	1.8560
8.0	1.0921	1.1333	1.1756	1.2189	1.2631	1.3082	1.3541	1.4009	1.4485	1.4968	1.5457	1.5953	1.6456	1.6964	1.7478	1.7997	1.8522
9.0	1.0906	1.1317	1.1738	1.2169	1.2610	1.3060	1.3518	1.3984	1.4458	1.4939	1.5427	1.5922	1.6423	1.6930	1.7442	1.7960	1.8483
10.0	1.0891	1.1301	1.1721	1.2150	1.2589	1.3037	1.3494	1.3959	1.4431	1.4911	1.5398	1.5891	1.6390	1.6895	1.7406	1.7923	1.8445
11.0	1.0876	1.1285	1.1703	1.2131	1.2569	1.3015	1.3471	1.3934	1.4405	1.4883	1.5368	1.5860	1.6358	1.6861	1.7371	1.7886	1.8407
12.0	1.0861	1.1268	1.1686	1.2112	1.2549	1.2994	1.3447	1.3909	1.4378	1.4855	1.5339	1.5829	1.6325	1.6828	1.7336	1.7849	1.8369
13.0	1.0846	1.1252	1.1668	1.2093	1.2528	1.2972	1.3424	1.3884	1.4352	1.4827	1.5309	1.5798	1.6293	1.6794	1.7300	1.7813	1.8331
14.0	1.0832	1.1236	1.1651	1.2075	1.2508	1.2950	1.3401	1.3859	1.4326	1.4799	1.5280	1.5767	1.6261	1.6760	1.7265	1.7776	1.8293
15.0	1.0817	1.1220	1.1633	1.2056	1.2488	1.2928	1.3377	1.3835	1.4300	1.4772	1.5251	1.5737	1.6229	1.6727	1.7231	1.7740	1.8256
16.0	1.0802	1.1204	1.1616	1.2037	1.2467	1.2907	1.3354	1.3810	1.4274	1.4744	1.5222	1.5706	1.6197	1.6693	1.7195	1.7704	1.8218
17.0	1.0788	1.1188	1.1599	1.2018	1.2447	1.2885	1.3331	1.3786	1.4248	1.4717	1.5193	1.5676	1.6165	1.6660	1.7161	1.7668	1.8181

Calculate diffusion ✕

Oxygen diffusion coefficient

Temperature (°C)

Salinity (‰)

D(10⁻⁵ cm² s⁻¹)

Apply
Cancel

Sediment diffusivity

- How should diffusivity be calculated?

DS and Theta

Settings
Zones with different diffusion parameters: 2

Calculate
D0 coefficient free water(10⁻⁵ cm² s⁻¹)
O2 Table 2.104

Sediment diffusion coefficient equation
DS = Por*D
DS = Por²*D
DS = D/(1+3*(1-Por))

Start	End	Porosity	DS
-1500	0	1	2.1
0	1200	0.8	1.683E-05

Confirm Cancel Enter DS manually

Fine sand/mud with high porosity, biofilm

Compact mud, low porosity

All kinds of sediments

Now a demonstration of the software
ProfilingDemo.upro – Profile 6.1

Porosity measurement

- Core sediment slicing e.g. every 5 mm
- Weight loss after drying / specific density water

Result – Statistics

Statistics

Save solution Export selected analysis

No. of Zones	SSE	P-Value	Top Conc (μmol/L)	Bottom Conc (μmol/L)	Top Flux (nmol cm ⁻² s ⁻¹)	Bottom Flux (nmol cm ⁻² s ⁻¹)	Integrated pr (nmol cm ⁻² s ⁻¹)
2	891.00	0.000	286.74	1.00	-0.015	0.000	0.015
3	688.98	0.071	284.07	1.00	-0.018	0.000	0.018
4	219.23	0.002	287.24	1.00	-0.014	0.000	0.014
5	102.52	0.007	282.56	1.00	-0.023	0.000	0.023
6	93.61	0.361	283.10	1.00	-0.021	0.000	0.021

- Choose the right number of zones:
 - P: Did adding this zone improve the fit ($P < 0.05$)
 - SSE: Sum of squared errors
 - SSE and P-value should be as small as possible.
 - Look at the results and use common sense


Result – Avoid obvious errors

Statistics


Save solution Export selected analysis

No. of Zones	SSE	P-Value	Top Conc (μmol/L)	Bottom Conc (μmol/L)	Top Flux (nmol cm ⁻² s ⁻¹)	Bottom Flux (nmol cm ⁻² s ⁻¹)	Integrated pr (nmol cm ⁻² s ⁻¹)
2	891.00	0.000	286.74	1.00	-0.015	0.000	0.015
3	688.98	0.071	284.07	1.00	-0.018	0.000	0.018
4	219.23	0.002	287.24	1.00	-0.014	0.000	0.014
5	102.52	0.007	282.56	1.00	-0.023	0.000	0.023
6	93.61	0.361	283.10	1.00	-0.021	0.000	0.021

- Production or consumption in zones where it is not possible
 - Analysis too deep into the anoxic zone
 - Maximum number of zones too high




Results – Number of zones



Statistics

Save solution Export selected analysis

No. of Zones	SSE	P-Value	Top Conc (µmol/L)	Bottom Conc (µmol/L)	Top Flux (nmol cm ⁻² s ⁻¹)	Bottom Flux (nmol cm ⁻² s ⁻¹)	Integrated pr (nmol cm ⁻² s ⁻¹)
6	7.58	0.069	165.46	0.27	0.058	0.000	-0.058
7	6.07	0.347	165.50	0.27	0.058	0.000	-0.058
8	6.16	1.000	165.51	0.27	0.058	0.000	-0.058
9	5.10	0.457	166.02	0.27	0.064	0.000	-0.064
10	6.00	1.000	165.99	0.27	0.063	0.000	-0.063



Statistics

Save solution Export selected analysis

No. of Zones	SSE	P-Value	Top Conc (µmol/L)	Bottom Conc (µmol/L)	Top Flux (nmol cm ⁻² s ⁻¹)	Bottom Flux (nmol cm ⁻² s ⁻¹)	Integrated pr (nmol cm ⁻² s ⁻¹)
3	37.81	0.212	167.63	0.27	0.066	0.000	-0.066
4	20.13	0.030	165.59	0.27	0.059	0.000	-0.059
5	13.34	0.102	165.06	0.27	0.054	0.000	-0.054
6	7.58	0.069	165.46	0.27	0.058	0.000	-0.058
7	6.07	0.347	165.50	0.27	0.058	0.000	-0.058



Other solutes than O₂

Boundary Conditions

Bottom conc + bottom flux


Top conc + bottom conc

Top conc + top flux

Top conc + bottom flux

Bottom conc + top flux


Bottom conc + bottom flux



Boundary conditions

Bottom conc. = 0

Bottom flux = 0



Boundary conditions

Top conc. = 0

Top flux = 0

Play around

- Boundary conditions
- Depth interval
- The oxygen diffusion rate - e.g. by changing the temperature and salinity
- Change the formula for the D_s calculation
- Porosity

Export data

Statistics

Save solution Export selected analysis

No. of Zones	SSE	P-Value	Top Conc (µmol/L)	Bottom Conc (µmol/L)	Top Flux (nmol cm ⁻² s ⁻¹)	Bottom Flux (nmol cm ⁻² s ⁻¹)	Integrated prod (nmol cm ⁻² s ⁻¹)
1	2018.07	0.000	147.95	0.27	0.044	0.000	-0.044
2	46.17	0.000	168.77	0.27	0.068	0.000	-0.068
3	15.29	0.007	165.90	0.27	0.061	0.000	-0.061
4	8.61	0.057	165.24	0.27	0.058	0.000	-0.058
5	7.68	0.444	165.70	0.27	0.061	0.000	-0.061

Choose how to export data

Export as Excel (.xlsx) file (Recommended).

Export as CSV files.

Ok Cancel **Note for csv: ',' is used as delimiter and '.' as decimal separator**

	A	B	C	D	E	F	G	H	I
1	Zone	SSE	P-Value	Top Concentration (umol/L)	Bottom Concentration (umol/L)	Top Flux (nmol cm ⁻² s ⁻¹)	Bottom Flux (nmol cm ⁻² s ⁻¹)	Integrated prod (nmol cm ⁻² s ⁻¹)	
2	1	2018.069859	0	147.9542313	0.270799994	0.044186883	-9.22382E-14	-0.044186883	
3	2	46.1681324	1.18735E-05	168.7682868	0.270799994	0.068284913	-1.43317E-13	-0.068284913	
4	3	15.28808306	0.006842208	165.9006005	0.270799994	0.060844866	-1.27845E-13	-0.060844866	
5	4	8.606835457	0.056830583	165.2440117	0.270799994	0.057699309	-1.1953E-13	-0.057699309	
6	5	7.684225285	0.444182386	165.6985895	0.270799994	0.061102369	-1.26345E-13	-0.061102369	
7									
8									
9									
10									
11									
12									
13									
14									
15									
16									
17									

Summary



SensorTrace Profiling – Activity calculations

- Calculate fluxes, consumption and production rates from high resolution concentration profiles.
- Fick's first and second law.
- Diffusivity and boundary conditions must be defined
- Use a stepwise optimization using the model-line and statistical values



Time for questions !

Unisense Microsensor Academy:
<https://www.unisense.com/support/knowledge>

Contact us: sales@unisense.com