



ENABLING
MICROSCALE
RESEARCH



Mitigation of N₂O emissions from wastewater biofilms

Microsensors confirm that counter-diffusion biofilms have lower N₂O emissions than co-diffusion biofilms



Introduction

N₂O is a very potent greenhouse gas (GHG) and accounts for up to 90% of the GHG emissions from wastewater treatment plants (WWTPs). N₂O is an intermediate product in biological treatment processes at WWTPs.

In this study, Professor Akihiko Terada and his research group at Tokyo University of Agriculture and Technology have investigated mitigation of N₂O emissions in a membrane-aerated biofilm reactor (MABR).

In a conventional biofilm reactor (CBR), the oxygen and electron donors (organic carbon and NH₄⁺) are supplied from the top of the biofilm from the liquid phase (co-diffusion). In an MABR, oxygen is supplied from the bottom of the biofilm through a gas-permeable membrane whereas the electron donors are supplied from the top of the biofilm (counter-diffusion). With this geometry, there will be a part in the middle of the MABR biofilm where electron acceptors co-exist with an electron donor, and this allows for simultaneous nitrification/denitrification, which could facilitate N₂O mitigation.

Laboratory setup

The Unisense MicroProfiling System was used to complete high resolution concentration profiles throughout the depth of the biofilms in co-diffusion and counter-diffusion biofilm reactors (Figure 1).

The biofilms were approximately 1500 μm thick. The researchers used an N₂O microsensor with a tip diameter of 25 μm (N₂O-25) and an O₂ microsensor with a tip diameter of 50 μm (OX-50) to make depth profiles throughout the biofilm inside of the biofilm reactor.




Figure 1: MicroProfiling setup showing the microsensor inserted through a port into the biofilm reactor.

Results and conclusion

The oxygen microprofiles in the counter-diffusion biofilm showed an oxygen penetration depth of 400 μm into the biofilm from the bottom and the O₂ concentration was highest at the biofilm-membrane interface (0 μm) where the air is supplied (Figure 2).

The N₂O concentration decreased just after O₂ depletion. The N₂O concentration at the biofilm-liquid interface was approximately 130 times lower in the MABR compared to the CBR.

From the concentration profiles, using the Fick's second law of diffusion, the researchers could calculate the N₂O production/consumption rates at the different depths in the biofilms (data not shown here). The authors found adjacent N₂O production/consumption hot spots and the positions of these most likely explained the increased N₂O consumption in the MABR biofilm.

The researchers could conclude that there was far less N₂O emission from the MABR compared to the conventional CBR and that the MABR is a promising technology for mitigation of N₂O emissions from WWTPs. You can read more in the article by Kinh et al. "Counter-diffusion biofilms have lower N₂O emissions than co-diffusion biofilms during simultaneous nitrification and denitrification: Insights from depth-profile analysis", Water Research 124 (2017) 363-371.




Figure 2: O₂ and N₂O concentration profiles within the MABR biofilm on day 95. The data are after Kinh et al. (2017).

Prof. Akihiko Terada says:

“Unisense O₂ and N₂O microsensors allow fast, accurate, and reliable activity measurements of microorganisms in suspensions and biofilms. They have provided our research group with opportunities to lead to exciting discoveries of bacteria and biofilm hotspots responsible for N₂O consumption.

Staffs are always kind and listen to our requests to improve/retrofit their products. We are sure to enjoy a scientific journey with Unisense microsensors as buddies”.

Suggested products



N₂O-25



MicroProfiling System



SensorTrace Profiling

Related publications

Qi et al. "Organic carbon determines nitrous oxide consumption activity of clade I and II nosZ bacteria: Genomic and biokinetic insights" Water Research Vol 209, 117910, 2022.

Suenaga et al. "Combination of 15N Tracer and Microbial Analyses Discloses N₂O Sink Potential of the Anammox Community", Environ. Sci. Technol. 2021, 55, 9231-9242

Suenaga et al. "Enrichment, Isolation, and Characterization of High-Affinity N₂O Reducing Bacteria in a Gas-Permeable Membrane Reactor" Environ. Sci. Technol. 2019

Suenaga et al. "Immobilization of Azospira sp. strain I13 by gel entrapment for mitigation of N₂O from biological wastewater treatment plants: Biokinetic characterization and modeling" Journal of Bioscience and Bioengineering Vol 126, No. 2, 213-219, 2018.

Jiang et al. "New insight into CO₂ -mediated denitrification process in H₂ -based membrane biofilm reactor: An experimental and modeling study", Water Research 184, 2020.

Lackner et al. Nitrification performance in membrane-aerated biofilm reactors differs from conventional biofilm systems, Water Research 44 (2010) 6073-6084