



GETTING STARTED GUIDE AND UNISENSE STANDARD SPECIFICATIONS FOR OXYGEN OPTODES (OPTO-SERIES)

Important! Test sensor upon receipt!

This sensor has been successfully tested prior to shipping, however some sensors suffer from rough transportation. Therefore, it is important that you test the sensor upon arrival.

REPLACEMENT OF DEFECTIVE SENSORS

Unisense will replace the sensor if it does not meet the specifications below, provided that:

1. A test is performed upon receipt without breaking the seal (Note! No seal on MR-sensors for testing purposes)
2. The complaint is given to Unisense **within two weeks** from receipt of the equipment.

GUARANTEED LIFETIME

Unisense guarantees the optode sensor a minimum lifetime of 1 year or a 1 million datapoints, whichever comes first, on condition of correct storage and use according to the manual.

INDIVIDUAL SENSOR CALIBRATION IS REQUIRED

Our sensors are handmade and as the sensor signal relies on the exact geometry of the sensor tip (micrometer scale), some variation must be expected.

SIGNAL AMPLIFICATION

Unisense oxygen optodes should be connected to a Opto-series or UniAmp-series amplifier.

STANDARD OXYGEN OPTODES ARE FUNCTIONING CORRECTLY IF (AT ROOM TEMPERATURE):

- The 90 % response time is within the specified range (see table below)
- The uncalibrated phase is within the specified range (see table below)

	OPTO-430	OPTO-430 FAST	OPTO-MR	OPTO-50	OPTO-3000
90 % response time (in sec.)	< 3 sec	< 0.3 sec	< 3 sec	< 5 sec	< 15 sec
100 % air saturation - phase	< 22	< 22	< 22	< 22	< 22
0 % oxygen - phase	> 46	> 46	> 46	> 40	> 46
Minimum datapoints (millions)	1.5	1	1.5	1	3



GETTING STARTED WITH OXYGEN OPTODES (OPTO-SERIES)

1. UNPACKING

- Remove the grey shock-absorbing plastic net and inspect the sensor visually. Leave the sensor in the protection tube for testing.

2. CONNECT THE OPTODE TO THE OPTO- OR UNIAMP-SERIES AMPLIFIER

- Remove the cap from the MicroOptode connector and from the amplifier connector.
- Insert the optode plug into the connector on the amplifier and turn gently clockwise until the plug is locked firmly.
- Insert the E²PROM connector on the optode into the E²PROM connector on the amplifier (no E²PROM on OPTO-Field sensors).



Unisense
MicroOptode

3. CALIBRATE THE SENSOR

- Use air saturated water as one calibration point. This is easily done using the CAL300 calibration chamber. The sensor may be dipped directly into the calibration chamber or the air saturated water may be injected into the protection tube using the calibration cap (see picture).
- Use the Unisense zero O₂ solution as the second calibration point. Inject the zero O₂ solution into the calibration cap and wait for the sensor to respond.
- For alternative calibration method, see the Oxygen MicroOptode - Opto Series User Manual.



CAL300 with
microsensors
and bubbling
with air.



Injecting calibration
liquid into protection
tube using the cali-
bration cap.

4. APPROVE THE SENSOR

- Compare the calibration values to the specifications given on the previous page. If necessary, see Troubleshooting in the optode manual or contact support (see below)

5. ADJUSTING THE MEASUREMENT FREQUENCY

- The lifespan of the optodes is determined by the amount of measurements. In order to prolong the lifespan of the optode, adjust the measuring frequency as required for your application.
- The frequency can be adjusted in the SensorTrace software

6. STORAGE

- When not in use unplug the optode and store it with the protection tube mounted at room temperature (preferable in low light or darkness).

USEFUL TOOLS



For support go to
www.unisense.com/support or
contact sales@unisense.com



Get the full manuals for all
sensors, equipment & software
at www.unisense.com/manuals



O₂ MicroOptode
Manual



Calkit-O₂
Manual



SensorTrace Suite
Manual



Find SDS for Cali-
bration Kit here