



ENABLING  
MICROSCALE  
RESEARCH



# Mitigation of N<sub>2</sub>O emissions from wastewater biofilms

Microsensors confirm that counter-diffusion biofilms have lower N<sub>2</sub>O emissions than co-diffusion biofilms



## Introduction

N<sub>2</sub>O is a very potent greenhouse gas (GHG) and accounts for up to 90% of the GHG emissions from wastewater treatment plants (WWTPs). N<sub>2</sub>O is an intermediate product in biological treatment processes at WWTPs.


In this study, Professor Akihiko Terada and his research group at Tokyo University of Agriculture and Technology have investigated mitigation of N<sub>2</sub>O emissions in a membrane-aerated biofilm reactor (MABR).

In a conventional biofilm reactor (CBR), the oxygen and electron donors (organic carbon and NH<sub>4</sub><sup>+</sup>) are supplied from the top of the biofilm from the liquid phase (co-diffusion). In an MABR, oxygen is supplied from the bottom of the biofilm through a gas-permeable membrane whereas the electron donors are supplied from the top of the biofilm (counter-diffusion). With this geometry, there will be a part in the middle of the MABR biofilm where electron acceptors co-exist with an electron donor. This allows for simultaneous nitrification/denitrification, which could facilitate N<sub>2</sub>O mitigation.

## Laboratory setup

The Unisense MicroProfiling System was used to complete high resolution concentration profiles throughout the depth of the biofilms in co-diffusion and counter-diffusion biofilm reactors (Figure 1).

The biofilms were approximately 1,500 μm thick. The researchers used an N<sub>2</sub>O microsensor with a tip diameter of 25 μm (N<sub>2</sub>O-25) and an O<sub>2</sub> microsensor with a tip diameter of 50 μm (OX-50) to make depth profiles throughout the biofilm inside of the biofilm reactor.



**Figure 1:** MicroProfiling setup showing the microsensor inserted through a port into the biofilm reactor.

## Results and conclusion


The oxygen microprofiles in the counter-diffusion biofilm showed an oxygen penetration depth of 400 μm into the biofilm from the bottom and the O<sub>2</sub> concentration was highest at the biofilm-membrane interface (0 μm) where the air is supplied (Figure 2).

The N<sub>2</sub>O concentration decreased just after O<sub>2</sub> depletion. The N<sub>2</sub>O concentration at the biofilm-liquid interface was approximately 130 times lower in the MABR compared to the CBR.

From the concentration profiles, using the Fick's second law of diffusion, the researchers could calculate the N<sub>2</sub>O production/consumption rates at the different depths in the biofilms (data not shown here). The authors found adjacent N<sub>2</sub>O production/consumption hot spots and the positions of these most likely explained the increased N<sub>2</sub>O consumption in the MABR biofilm.

The researchers could conclude that there was far less N<sub>2</sub>O emission from the MABR compared to the conventional CBR and that the MABR is a promising technology for mitigation of N<sub>2</sub>O emissions from WWTPs.

You can read more in the article by Kinh et al. "Counter-diffusion biofilms have lower N<sub>2</sub>O emissions than co-diffusion biofilms during simultaneous nitrification and denitrification: Insights from depth-profile analysis", *Water Research* 124 (2017) 363-371.



**Figure 2:** O<sub>2</sub> and N<sub>2</sub>O concentration profiles within the MABR biofilm on day 95. The data are after Kinh et al. (2017).

**Prof. Akihiko Terada says:**

"Unisense O<sub>2</sub> and N<sub>2</sub>O microsensors allow fast, accurate, and reliable activity measurements of microorganisms in suspensions and biofilms. They have provided our research group with opportunities to lead to exciting discoveries of bacteria and biofilm hotspots responsible for N<sub>2</sub>O consumption.

Staff is always kind and listen to our requests to improve/retrofit their products. We are sure to enjoy a scientific journey with Unisense microsensors as buddies".

**Suggested products**



N<sub>2</sub>O-25



MicroProfiling System



SensorTrace Profiling

**Related publications**

Qi et al. "Organic carbon determines nitrous oxide consumption activity of clade I and II nosZ bacteria: Genomic and biokinetic insights" *Water Research* Vol 209, 117910, 2022.

Suenaga et al. "Combination of 15N Tracer and Microbial Analyses Discloses N<sub>2</sub>O Sink Potential of the Anammox Community", *Environ. Sci. Technol.* 2021, 55, 9231-9242

Suenaga et al. "Enrichment, Isolation, and Characterization of High-Affinity N<sub>2</sub>O Reducing Bacteria in a Gas-Permeable Membrane Reactor" *Environ. Sci. Technol.* 2019

Suenaga et al. "Immobilization of Azospira sp. strain I13 by gel entrapment for mitigation of N<sub>2</sub>O from biological wastewater treatment plants: Biokinetic characterization and modeling" *Journal of Bioscience and Bioengineering* Vol 126, No. 2, 213-219, 2018.

Jiang et al. "New insight into CO<sub>2</sub> -mediated denitrification process in H<sub>2</sub> -based membrane biofilm reactor: An experimental and modeling study", *Water Research* 184, 2020.

Lackner et al. Nitration performance in membrane-aerated biofilm reactors differs from conventional biofilm systems, *Water Research* 44 (2010) 6073-6084